Ω	Ω CALIBRES NO. 35		CS08-19-C-035-E	4
OMEGA	Made by: pelrom	Date: 13.10.2004		

Replacement of movement calibre 1459 with calibre 1459B

PROBLEM

Calibre 1459 and its components can only be delivered while stocks last.

SOLUTION

The calibre 1459 movement is now being replaced by the calibre 1459B. The external dimensions and those of the casing are identical. Although the height of the movement is different when no battery is fitted, only the height with a battery fitted (1.40 mm) is of importance, and this is identical.

By using the same 9951 battery, the battery life is increased from 27 to 34 months.

CHARACTERISTICS

The calibre 1459 is fully exchangeable with the calibre 1459B. For details of the components that have been modified, please consult the 1459B technical manual.

<u>PS</u>: Movements to be exchanged must be returned to Omega SA, Spare Parts International, CH-2500 Bienne 4, Switzerland.